

# Anti-H<sub>8</sub> ((A/pintail duck/Alberta/114/1979(H<sub>8</sub>N<sub>4</sub>)))



---

Product Name	Anti-H <sub>8</sub> (A/pintail duck/Alberta/114/1979(H <sub>8</sub> N <sub>4</sub> ))
Catalog No	IT-003-033
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	H <sub>8</sub> (A/pintail duck/Alberta/114/1979(H <sub>8</sub> N <sub>4</sub> )) protein (Genebank accession <a href="#">ABB87729</a> )
Applications	Western Blotting (1:500-1:5000), ELISA(1:100-1:100000)
Specificity	Reacts with H <sub>8</sub> (H <sub>8</sub> N <sub>4</sub> ) protein. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1 (A/pintail duck/Alberta/114/1979(H<sub>8</sub>N<sub>4</sub>))