



Anti-HA(B/Victoria/2/87)

Product Name	Anti- HA(B/Victoria/2/87)
Catalog No	IT-003-032M5
Source	Mouse monoclonal antibody, clone 3C5
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant HA(Δ TM)(B/Victoria/2/87)(IT-001-B6p) protein (a.a. 16-548)(Genebank No. AAA43697).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HA(B/Victoria/2/87) only as tested so far. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-032M5

