

Anti-N7(H7N7)



Product Name	Anti-N2 (H7N7)
Catalog No	IT-003-028M2
Source	Mouse monoclonal antibody, clone 43E4
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant NA protein (A/chicken/Netherlands/1/03(H7N7)) (Genebank accession No. AAR11365)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most contemporary N7 (H7N7). Cross-reacts with N2 (H3N2). No cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

