

Anti-N2(H9N2)



Product Name	Anti-N2 (H9N2)
Catalog No	IT-003-027M2
Source	Mouse monoclonal antibody, clone 173A10
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant NA protein (A/chicken/HongKong/NT366/03(H9N2)) (Genebank accession No. AAV28392)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with most contemporary N2 (H9N2). Cross-reacts with N2 (H3N2). Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

