

Anti-N2(H9N2)



Product Name	Anti-N2 (H9N2)
Catalog No	IT-003-027M1
Source	Mouse monoclonal antibody, clone 187G7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant NA protein (A/chicken/HongKong/NT366/03(H9N2)) (Genebank accession No. AAV28392)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most contemporary N2 (H9N2). No cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

IT-003-027M1

