

# Anti-M1(A/California/04/2009(H1N1))



Product Name	Anti- M1(A/California/04/2009(H1N1))
Catalog No	IT-003-036M1
Source	Mouse monoclonal antibody, clone 7G2
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing M1(A/California/04/2009(H1N1)) (Genebank No. ACP44152)
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with M1 (A/California/04/2009(H1N1)) and related M1 from H1N1 subtypes. Cross-react with other M1 from other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

