


Anti-NS1 (A/California/06/2009/(H₁N₁))



Product Name	Anti-NS1 ((A/California/06/2009/(H ₁ N ₁)))
Catalog No	IT-003-037
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	NS1 ((A/California/06/2009/(H ₁ N ₁))) protein (Genebank accession ACP52562)
Applications	Western Blotting (1:200-1:1000), ELISA(1:100-1:10000)
Specificity	Reacts with NS1(H ₁ N ₁) protein. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1 ((A/California/06/2009/(H₁N₁)))