


# Anti-NA (New Caledonia/20/1999 (H<sub>1</sub>N<sub>1</sub>))



---

Product Name	Anti-NA(New Caledonia/20/1999/(H <sub>1</sub> N <sub>1</sub> ))
Catalog No	IT-003-039
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	NA(New Caledonia/20/1999/(H <sub>1</sub> N <sub>1</sub> )) protein (Genebank accession ABQ10080)
Applications	Western Blotting (1:200-1:1000), ELISA(1:100-1:10000)
Specificity	Reacts with NA(H <sub>1</sub> N <sub>1</sub> ) protein. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NA(New Caledonia/20/1999/(H<sub>1</sub>N<sub>1</sub>))