

# Anti-H<sub>1</sub>(H<sub>1</sub>N<sub>1</sub>)(A/Brisbane/59/2007)



---

Product Name	Anti-H <sub>1</sub> (Bris/2007)
Catalog No	IT-003-Ref-H1(Bri/2007)
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H <sub>1</sub> (H <sub>1</sub> N <sub>1</sub> ) (A/Brisbane/59/2007) (Genebank Accession # ACA28846)
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H <sub>1</sub> (H <sub>1</sub> N <sub>1</sub> ) (A/Brisbane/59/2007) and related H <sub>1</sub> . Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.