

Anti-H₁(H₁N₁)(A/New Caledonia/20/99)



Product Name	Anti-H ₁ (NC/1999)
Catalog No	IT-003-Ref-H ₁ (NC/1999)
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H ₁ (H ₁ N ₁) (A/New Caledonia/20/99) (Genebank Accession # AAP34324)
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H ₁ (H ₁ N ₁) (A/New Caledonia/20/99) and related H ₁ . Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.