

Anti-H₁(H₁N₁)(A/South Carolina/1/18)



Product Name	Anti-H ₁ (SC/1918)
Catalog No	IT-003-Ref-H1(SC/1918)
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H ₁ (H ₁ N ₁) (A/ South Carolina/1/18) (Genebank Accession # AAD17229)
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H ₁ (H ₁ N ₁) (A/ South Carolina/1/18) and related H ₁ . Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.