

## Anti-H<sub>1</sub>(H<sub>1</sub>N<sub>1</sub>)(A/Ohio/Swine/01/2007)



---

Product Name	Anti-H <sub>1</sub> (SW/2007)
Catalog No	IT-003-Ref-H1(SW/2007)
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H <sub>1</sub> (H <sub>1</sub> N <sub>1</sub> ) (A/Ohio/Swine/01/2007) (Genebank Accession # <a href="#">ACQ73203</a> )
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H <sub>1</sub> (H <sub>1</sub> N <sub>1</sub> ) (A/Ohio/Swine/01/2007) and current circulating H <sub>1</sub> from “Novel 2009 H <sub>1</sub> N <sub>1</sub> ”. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.