

Anti-H₁₂(H₁₂N₅)(A/green-winged teal/ALB/199/1991)



Product Name	Anti-H ₁₂ (H ₁₂ N ₅)(A/green-winged teal/ALB/199/1991)
Catalog No	IT-003-Ref-H ₁₂
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H ₁₂ (H ₁₂ N ₅)(A/green-winged teal/ALB/199/1991) (Genbank Accession # ABB88110)
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H ₁₂ (H ₁₂ N ₅) (A/green-winged teal/ALB/199/1991) and related H ₁₂ . Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.