

Anti-H₁₃(H₁₃N₈)(A/black-headed gull/Netherlands/1/00)



Product Name	Anti-H ₁₃ (H ₁₃ N ₈) (A/black-headed gull/Netherlands/1/00)
Catalog No	IT-003-Ref-H13
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H ₁₃ (H ₁₃ N ₈) (A/black-headed gull/Netherlands/1/00) (Genebank Accession # AAV91212)
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H ₁₃ (H ₁₃ N ₈) (A/black-headed gull/Netherlands/1/00) and related H ₁₃ viruses. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.