

Anti-H₃(H₃N₂)(A/Wisconsin/67/X-161/2005)



Product Name	Anti-H3(Wis/67/2005)
Catalog No	IT-003-Ref-H3(Wis/2005)
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H ₃ (H ₃ N ₂)(A/ Wisconsin/67/X-161/2005) (Genebank Accession # ABO37609)
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H ₃ (H ₃ N ₂)(A/ Wisconsin/67/X-161/2005) and related H ₃ . Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.