

# Anti-H<sub>3</sub>(H<sub>3</sub>N<sub>2</sub>)( A/Wyoming/3/03)



---

Product Name	Anti-H3(Wyoming/2003)
Catalog No	IT-003-Ref-H3(Wyo/2003)
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS
Immunogen	DNA vaccine expressing H <sub>3</sub> (H <sub>3</sub> N <sub>2</sub> )( A/Wyoming/3/03) (Genebank Accession # <a href="#">AAT08000</a> )
Applications	Neutralization assay; ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with H <sub>3</sub> (H <sub>3</sub> N <sub>2</sub> )( A/Wyoming/3/03) and related H <sub>3</sub> . Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.