

Anti-H1(A/swine/England/WVL7/1992(H1N1))



Product Name	Anti- H1(A/swine/England/WVL7/1992(H1N1))
Catalog No	IT-003-SW8M2
Source	Mouse monoclonal antibody, clone 350D9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing H1(A/swine/England/WVL7/1992(H1N1)) (Genebank No. ACO25133)
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H1(A/swine/England/WVL7/1992(H1N1)), H1(A/South Carolina/1/18/(H1N1)), H1(A/swine/Guangxi/13/ 2006(H1N2)) (weak). Not cross-react with other H1 and other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-SW8M2

