



# Anti-H1(A/Swine/(H1N1))

Product Name	Anti- H1(A/Swine/(H1N1))
Catalog No	IT-003-SWM1
Source	Mouse monoclonal antibody, clone 15B10
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing H1(A/California/06/2009(H1N1)) (Genebank No. <a href="#">ACP41935</a> )
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H1(A/California/06/2009(H1N1)). Cross-reacts with H1 (A/swine/Guangxi/13/2006(H1N2) and H1(A/New Jersey/8/1976 (H1N1))(weak). Not cross-react with other seasonal H1 and other subtypes.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-003-SWM1

