



Anti-gG (HSV-2)

Product Name	Anti-gG (HSV-2)
Catalog No	IT-005-010
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Human HSV-2 gG protein (a.a. 23-650) (Genebank accession #NP_044534)).
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HSV-2 gG proteins. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

