

Anti-gH/gL/gp42 Complex(EBV)



Product Name	Anti-gH/gL/gp42 Complex(EBV)
Catalog No	IT-005-033
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing ectodomains of gH/gL/gp42 complex (gH: a.a. 1-344, gL: a.a 1-136, gp42: a.a. 1-217) proteins (strain B95-8) (Genebank No. YP_401700)
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with gH/gL/gp42 complex proteins. Minimal cross-reactivity to other herpesviruses.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-005-033

