




Anti-E (West Nile Virus)

Product Name	Anti-E (WNV)
Catalog No	IT-006-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µl rabbit anti-sera with 0.1% sodium Azide.
Immunogen	E protein (West Nile virus (lineage I strain NY99)) (Genebank No. NC_009942)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with E protein or its subdomain (DIII) of West Nile Virus.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
1: 293 cell expressing E(DIII)
2: 293 cell expressing E ectodomain

Note: