



Anti-PreM/M (West Nile Virus)

Product Name	Anti-PreM/M (WNV)
Catalog No	IT-006-002
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µl of rabbit anti-sera with 0.1% sodium Azide.
Immunogen	PreM protein (West Nile virus (lineage I strain NY99)) (Genebank No. NC_009942)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with PreM/M protein of West Nile Virus.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing PreM protein

Note: