



# Anti-NS1(West Nile Virus)

---

Product Name	Anti-NS1 (WNV)
Catalog No	IT-006-003
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µl of rabbit sera with 0.1% sodium Azide.
Immunogen	NS1 protein (West Nile virus (lineage I strain NY99)) (Genebank No. NC_009942)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with NS1 protein of West Nile Virus.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

