

Anti-NS1(Dengue Virus 3)



Product Name	Anti-NS1(DENV3)
Catalog No	IT-006-0052M1
Source	Mouse monoclonal antibody, clone 5D4
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Dengue virus 3/KPS-4-0657/207) (a.a. 1-373) (Genebank No. AAW84288)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Cross-reacts with NS1 from DENV1, 2, JEV, WNV.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-006-0052M1

