

Anti-NS1(WNV)



Product Name	Anti- NS1(WNV)
Catalog No	IT-006-0053EM1
Source	Mouse monoclonal antibody, clone 6G7
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant NS1 (WNV) (strain HNY1999) (a.a. 782-1102) (Genebank Accession No. AAF18443).(IT-006-0053Ep)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with NS1 (WNV). Cross-reactivity not determined.
Storage	Stable for 6 months at 4°C from the date of shipment. For long term storage, aliquot and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

