




Anti-NP(Measles Virus)

| | |
|----------------|--|
| Product Name | Anti-NP(Measles virus) |
| Catalog No | IT-007-010 |
| Source | Rabbit Polyclonal Antibody |
| Formulation | Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide. |
| Immunogen | NP protein (Measles Virus)(Genebank No. ABV69625) |
| Applications | Western Blotting (1:200-1:1000) |
| Specificity | Reacts with NP protein from Measles Virus. Cross-reactivity not tested. |
| Storage | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



C: 293 cell extract control
E: 293 cell expressing NP(Measles virus)

Note: