



## E (DIII)(West Nile Virus)

---

Product Name	E (Domain III)(WNV)
Catalog No	IT-006-0011p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged West Nile Virus E protein ectodomain (a.a.590-696) (Genebank No. AAF18443).
Applications	WB standard, antibody ELISA, vaccine, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: E (DIII)(WNV) protein

Note: