

NS1 (DV4)



Product Name	NS1 (DV4)
Catalog No	IT-006-0055p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Dengue Virus-4 NS1 protein + V5 (Genebank No. ACH61692).
Applications	WB standard, antibody ELISA, vaccine, etc
Purity	>90%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1 (DV4) protein

Note: