



gG (HSV-1)

Product Name	gG (HSV-1)
Catalog No	IT-005-002p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged gG (HSV-1) protein (a.a. 34-200) (Genebank accession #CAF24800). A GCN4 trimeric domain was fused to its C-terminus. The protein predominantly exists as trimer.
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified gG (HSV-1)/GCN4 fusion protein.