

gC (Ectodomain)(HSV-2)



Product Name	gC(Ectodomain)(HSV-2)
Catalog No	IT-005-011p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged gC(Ectodomain)(HSV-2)(a.a.28-445) (Genebank #: AFM93864).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

