


Human Complement Receptor Type 2 (CR2/CD21) Ectodomain



Product Name	Human Complement Receptor Type 2 (CR2/CD21) Ectodomain
Catalog No	IT-005-030CR2(Δ TM)
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged Human Complement Receptor Type 2 (CR2/CD21) Ectodomain(a.a.21-1035)(Genebank Accession No. BAD97347). hCR2 is the main receptor for EBV infection of human B cells.
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Complement Receptor Type 2 (CR2/CD21) Ectodomain