

# gD (Ectodomain)(HSV-1)



---

Product Name	gD(Ectodomain)(HSV-1)
Catalog No	IT-005-055p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged gD(Ectodomain)(HSV-1) (a.a.27-310) (Genebank #: AAK19597).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

