


NS3 (HCV/Genotype-3a/Con)



Product Name	NS3 (HCV/Genotype-3a/Con)
Catalog No	IT-004-071-3aEp
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS with 8M urea.
Protein	GST-fused HCV NS3 Genotype-3a consensus immunodominant regions (aa 1199-1466).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: GST-NS3 (HCV/Genotype-3a/Con)

Protein SEQ:

alqfipvetlstqarspsfsdnstppavpqsyqvgylhaptgsgkstkvpaaayvaqgynvlvlnpsvaatlgfgsfmsraygi
dpnirtgnrtvtgaktystygfkladggcsggaydviiicdechaqdatsilgigtvldqaetagvrltvlatappgsitvphsn
ieevalgsegeipfygkaipiallkggrhlifchskkkcdeiasklrgmglnavayyrgldvsviptgdvvvcatdalmtgft
gdfdsvidcnvave