

NS3 (HCV/Genotype-6)



Product Name	NS3 (HCV/Genotype-6)
Catalog No	IT-004-075Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in 1.5M urea, 25mM Tris-HCl pH-8, 0.2% Triton-X & 50% Glycerol. Reconstitute the protein with 100 µl of Millipore water.
Protein	The E.coli derived recombinant protein contains the HCV NS3 immunodominant regions, amino acids 1356-1459. The protein is fused to a GST tag at N-Terminus.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.