

# NS5 (HCV/Genotype-1)

---



Product Name	NS5 (HCV/Genotype-1)
Catalog No	IT-004-077Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in 50mM Tris. pH 8 & 5mM EDTA. Reconstitute the protein with 100 µl of Millipore water.
Protein	The E.coli derived recombinant protein contains the HCV NS5 Genotype1 immunodominant regions.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.