

gp140(HIV-1/Clade B'/C(CN54))



Product Name	gp140(HIV-1/Clade B'/C(CN54))
Catalog No	IT-001-161p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp140 (HIV-1/Clade B'/C(CN54)) protein (Genebank #: AX149771) (a.a. 29-676)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140 (HIV-1/Clade B'/C(CN54))