

gp140(HIV-1/Clade B'/C)(Consensus)



Product Name	gp140(HIV-1/Clade B'/C) (Consensus)
Catalog No	IT-001-162p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp140 (HIV-1/Clade B'/C)(Consensus) protein (http://www.hiv.lanl.gov/cIT-bin/NEWALIGN/align.cgi) (a.a. 35-674). It has 99% homology to a circulating strain (97CNGX-7F)(Genebank #: AAG38915) isolated in China in 2000.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140 (HIV-1/Clade B'/C)(Consensus)