



hCD317(Δ TMs)

Product Name	hCD317(Δ TMs)
Catalog No	IT-001-317p
Source	Purified protein from transfected 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged ectodomain of human CD317 (Tetherin) protein (a.a.24-162). CD317 has been recently demonstrated to cause retention of retroviral virions on cell surfaces (Nature 451 (7177), 425-430 (2008))
Applications	Interference of Vpr/CD317 interaction, WB standard, , etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

