



# Human CD4(D1D2)

---

Product Name	Human CD4(D1D2 region)
Catalog No	IT-001-CD4d1d2p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged human CD4(D1D2) region (a.a.26-206 of CD4 precursor) (Genebank No. NP_000607)
Applications	HIV-1 entry inhibitor, WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human CD4(D1D2 region)

Note: