



gp120(Con_A2)(HIV-1)

Product Name	gp120(Con_A2)
Catalog No	IT-001-CONA2p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (Con_A2/2006) protein (a.a.34-518) http://www.hiv.lanl.gov/cgi-bin/CONSENSUS_DOWNLOAD/ConsensusDownloader.cgi .
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

