




# gp120(pTHRO4156)(HIV-1)

---

|                |  |
|----------------|--|
| Product Name   | gp120(pTHRO4156)(Clade B)  |
| Catalog No     | IT-001-RB10p   |
| Source         | Viral protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.  |
| Protein        | 6xHis tagged HIV-1 gp120(pTHRO4156)(Clade B) protein (a.a.34-518) (Genebank No. AY835448).   |
| Applications   | WB standard, antibody ELISA, etc   |
| Purity         | >95%   |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |



E: gp120(pTHRO4156) protein

Note: