

HRP-Goat anti-Mouse IgG(H+L)



Product Name	HRP-Goat anti-Mouse IgG(H+L)
Catalog No	IT-200-010
Source	Antigen-specific-affinity purified IgG from goat anti-serum
Formulation	Each vial contains 0.5 ml of conjugated antibody (0.8 mg/ml) in PBS(pH7.6) with 15 mg/ml IgG-free BSA.
Immunogen	Purified mouse IgG
Applications	WB(1:10,000- 1:200,000), ELISA(1:5,000-1:100,000), immunohistochemistry(1:500-5,000)
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.