

HRP-Goat anti-Rabbit IgG(H+L)



Product Name	HRP-Goat anti-Rabbit IgG(H+L)
Catalog No	IT-200-011
Source	Antigen-specific-affinity purified IgG from goat anti-serum
Formulation	Each vial contains 0.5 ml of conjugated antibody (0.8 mg/ml) in PBS(pH7.6) with 15 mg/ml IgG-free BSA.
Immunogen	Purified rabbit IgG
Applications	WB(1:10,000- 1:200,000), ELISA(1:5,000-1:100,000), immunohistochemistry(1:500-5,000)
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.