


VP40 (Marburgvirus/Lake Victoria)



Product Name	VP40 (Marburgvirus/Lake Victoria)
Catalog No	IT-014-012Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged Matrix protein VP40 (Marburgvirus/Lake Victoria) (GenBank # Q1PD51, aa 1-303)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: VP40 (Marburgvirus/Lake Victoria)