


Spike Protein (MERS-CoV/EMC 2C/2012)



Product Name	Spike Protein (MERS-CoV/EMC 2C/2012)
Catalog No	IT-002-010p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Spike protein (MERS-Cov/EMC 2C/2012)(a.a.18-725) (Genebank No. AFS88936).
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Spike protein (MERS-CoV/EMC 2C/2012)

Note: