

Antigen TB15.3 (M. Tuberculosis/ CDC1551)



Product Name	Antigen TB15.3 (M. Tuberculosis/CDC1551)
Catalog No	IT-500-019Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged TB15.3/Rv1636/MT1672 protein (a.a. 1-146) (GenBank Accession No. Y1636_MYCTO).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

