

HA1(A/tufted duck/Germany-SH/R8444/2016) (H5N8)



Product Name	HA1(A/tufted duck/Germany-SH/R8444/2016)(H5N8)
Catalog No	IT-003-00536p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1 (A/tufted duck/Germany-SH/R8444/2016) (H5N8) (a.a.18-341) (GISAID No. EPI860509).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/tufted duck/Germany-SH/R8444/2016)(H5N8)

Note: