

HA(Δ TM)(A/Slovenia/2903/2015)(H1N1)



Product Name	HA(Δ TM)(A/Slovenia/2903/2015)(H1N1)
Catalog No	IT-003-00107 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA(A/Slovenia/2903/2015)(H1N1)(a.a.18-529) (GISAID No. EPI768541).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq 90\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

