




# Human GITRL

---

Product Name	Human GITRL
Catalog No	IT-300-031p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Human GITRL (a.a.72-199) (GenBank# NP_005083).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human GITRL

Note: