



Human CTLA-4 / CD152

Product Name	Human CTLA-4 / CD152
Catalog No	IT-300-032p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Human CTLA-4 (a.a.36-161) (GenBank# NP_005205).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

